


DQEPro

Radiography DQE Sample Report

- Wireless panel with poor performance due to short integration window
- Clinical-environment testing



This report was generated automatically by *DQEPro*. It summarizes the MTF and DQE performance of a clinical CsI-based flat-panel detector. The full DQE assessment was completed in 15 minutes in the x-ray room of a busy clinical department.

Observations:

1. The linearity test showed the system saturated at 75 μGy (8.6 mR), corresponding to a 1-sec exposure. This is a wireless detector with a preset integration window of 0.5 sec (a common value). The detector stops acquiring the image at 0.5 sec causing the pixel value to be flat (at 1290) even as the exposure continues. This is a serious problem with many wireless systems.

Conclusions:

Wireless detectors must be checked for adequate exposure window settings. In this example, the default window was too short, resulting in poor image quality and high patient exposures with large patients. This problem is particularly acute with low-mA generators. It can be avoided by increasing the exposure window setting.

X-ray spectrum:

Half-value layer: 7.1 mm Al

Best fit: SaturatedLinear

Nominal image air KERMA, exposure: 1.05 uGy, 0.120 mR

Average open pixel value (DV): 0.0

Average dark pixel value (DV): DVd = 0.0

Estimated saturation exposure time: 0.9 sec

Estimated saturation exposure: 75.0 uGy 8.6 mR


This system has a saturated linear response. Linearization is not possible.

Linearization accuracy will be improved by modifying the protocol to ensure at least three exposure levels are below the nominal exposure (0.12 mR, 1.0 uGy).


The measured slope (gradient) of the system response is neither all positive nor all negative. This generally means non-linear post-processing such as auto windowing or histogram equalization has been applied to the image.

Detector output saturates at approximately 0.9 sec of exposure (75.0 uGy, 8.55 mR). This may be due to inadequate detector integration window and will prevent correct operation with large patients. DQE cannot be determined. Increasing integration window to 3.3 sec or more may solve this problem.


Generator Linearity



System Response Measured



System Gradient Measured



X-ray spectrum:

Half-value layer: 7.1 mm Al

Number of ignored pixels in ROI signal calculation: 0 of 860160

Dark noise variance from analysis: 0.0

Dark noise variance from analysis if linearized: 0.0

Quantum-noise limit exposure: could not be determined

This system saturated (overexposure) and therefore the quantum-noise limit can not be determined.

Two or more images were acquired at only 0 different non-zero mAs settings in the linearity test. This is too few to complete the quantum noise analysis. Two or more images must be acquired at each of at least 5 different non-zero mAs settings.

The quantum-noise limit could not be determined because no dark images were acquired. To determine the quantum noise limit, ensure at least 2 dark images or 2 x-edge images are acquired.

The quantum limit analysis could not be completed because dark image pixel variance is zero.

The quantum-noise limit can not be determined because measured dark noise is greater than the lowest-mAs noise. This is not physical and could be caused by detector instability or some other problem with image data.

Quantum noise limit could not be determined. Be sure your images are free of artifacts and at least two images are acquired at each exposure level.

Data folder: \Users\DQE\data

These waveforms show x-ray exposure profiles measured by DQEPro. Green line segments correspond to identified exposure pulses. Leading and trailing baselines are dark blue and inter-pulse baselines are light blue. The first pulse should start at 0 sec. Ensure pulses are correctly identified and free of undesirable artifacts that may indicate problems with the x-ray generator or tube. Exposure values correspond to image plane based on Pro-image distance entered on panel. Double-click on any waveform to open a review window.




Image data folder: \Users\DQE\data\images

FileName	Format	Width	Height	Depth	bits	kV	mA	msec	mAs	SID	Pixel	um
1: 11_45_212.raw	vw2560x307	[2560	3072	1]	0	0	0	0	0.0	0	[0	0]
2: 11_45_213.raw	vw2560x307	[2560	3072	1]	0	0	0	0	0.0	0	[0	0]
3: 11_45_214.raw	vw2560x307	[2560	3072	1]	0	0	0	0	0.0	0	[0	0]
4: 11_45_215.raw	vw2560x307	[2560	3072	1]	0	0	0	0	0.0	0	[0	0]
5: 11_45_216.raw	vw2560x307	[2560	3072	1]	0	0	0	0	0.0	0	[0	0]
6: 11_45_217.raw	vw2560x307	[2560	3072	1]	0	0	0	0	0.0	0	[0	0]
7: 11_45_218.raw	vw2560x307	[2560	3072	1]	0	0	0	0	0.0	0	[0	0]
8: 11_45_219.raw	vw2560x307	[2560	3072	1]	0	0	0	0	0.0	0	[0	0]


Image data folder: \Users\DQE\data

These images are thumbnail illustrations of each image. For dynamic systems, only the first few images are shown for each exposure. Ensure each image is free of undesirable artifacts that may indicate problems with the x-ray system. Double-click on thumbnail to open review window.

1: 11_45_212.raw




2: 11_45_213.raw



3: 11_45_214.raw




4: 11_45_215.raw



5: 11_45_216.raw



6: 11_45_217.raw




7: 11_45_218.raw



8: 11_45_219.raw



X-ray spectrum: Radiography, RQA-5 (70 kV)
Target material: W
Set kV: 74
Set mAs: 1.25
Half-value layer: 7.1 mmAl
Mean photon energy: 53.1 keV
Spectrum Q_o: 30.2 q/mm²/nGy



Study comment: Sample study

DQEPro acquisition version: 4.3.4-rc1, Mac64 10.8 (Mountainlion) hostname username

DQEPro analysis version: 4.7.0-rc7, Windows-64bit Windows 10 Cunningham-spectre icunning

DQEPro hardware version: 2.3, 000000

DQEPro firmware version: 2.18.5

Grid:

X-ray spectrum: Radiography, RQA-5 (70 kV)

Half-value layer: 7.1 mmAl

Set technique: 74 kV, mA not specified, 1.3 mAs

Air Temperature and Pressure: 24.9 C, 103.8 kPa

Waveform sampling rate: 4.0 kHz

Image-plane exposure, air KERMA: 0.120 mR, 1.05 uGy

Source-image distance: 143.0 cm (DQEPro)

Specified DQEPro-image distance: 1.0 cm

DQEPro-image distance: 1.0 cm (Specified by user)

Pixel size from image headers, x y: NA

System response: SaturatedLinear

Data Folder: \Users\DQE\data

Image Folder: \Users\DQE\data\images

Warnings:

- 1: Linearization accuracy will be improved by modifying the protocol to ensure at least three exposure levels are below the nominal exposure (0.12 mR, 1.0 uGy).
- 2: The measured slope (gradient) of the system response is neither all positive nor all negative. This generally means non-linear post-processing such as auto windowing or histogram equalization has been applied to the image.
- 3: Detector output saturates at approximately 0.9 sec of exposure (75.0 uGy, 8.55 mR). This may be due to inadequate detector integration window and will prevent correct operation with large patients. DQE cannot be determined. Increasing integration window to 3.3 sec or more may solve this problem.
- 4: This system saturated (overexposure) and therefore the quantum-noise limit can not be determined.
- 5: Two or more images were acquired at only 0 different non-zero mAs settings in the linearity test. This is too few to complete the quantum noise analysis. Two or more images must be acquired at each of at least 5 different non-zero mAs settings.
- 6: The quantum-noise limit could not be determined because no dark images were acquired. To determine the quantum noise limit, ensure at least 2 dark images or 2 x-edge images are acquired.
- 7: The quantum limit analysis could not be completed because dark image pixel variance is zero.
- 8: The quantum-noise limit can not be determined because measured dark noise is greater than the lowest-mAs noise. This is not physical and could be caused by detector instability or some other problem with image data.
- 9: Quantum noise limit could not be determined. Be sure your images are free of artifacts and at least two images are acquired at each exposure level.
- 10: Images were saved in a non-dicom format. This means we cannot confirm "For Processing" or "Raw" images are used and that only linear post processing has been applied to image data according to NEMA standard XR 27-2012. The user must ensure only linear processing has been used to ensure IEC accuracy in results. In particular, noise-suppression algorithms and lossy compression may cause unpredictable results and must not be used.